## How to Create a Search Filter

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

The screenshot shows a 'Create search filter' dialog box with the following elements:


- Folder Information:**
  - Folder Name: [Text input field]
  - Description: [Text input field]
- Options:**
  - Link to OnCall
  - Auto Download Priors
  - Study DTTM ascending
- Search Conditions Table:**

*Study Date	Within last 30 days
Modality	*Any Modality
*Patient ID	*Any Patient ID
*Patient Name	*Any Patient Name
Specialty	*Any Specialty
- Show exam on top:**
  - ITR EMERGENCY STAT
  - High
  - Low
  - Medium
  - Routine
- Buttons:** Add condition, Delete condition, OK, Cancel

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Study DTTM ascending

\*Patient Name

- \*Patient Name
- \*Study Date
- \*Study key
- Acquisition date
- Approver
- Approver2
- Attending Doctor
- Body part
- Consulting Doctor
- Department(Branch of Service)
- Dictator
- Ed Status
- ER Doctor
- Exam Status**
- ICD-10 code
- Institution
- Modality
- Number of image count
- Operator
- Other patient ID
- Other patient name
- Passcode
- Patient Age
- Patient Birth Date
- Patient Location
- Patient Residency

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Routine

\*Any Exam Status

- \*Any Exam Status
- Examined
- Verified
- Dictating
- Unread
- Dictated
- Transcribing
- Transcribed
- Approved**
- Approved(2)
- Addendum\_dictating
- Addendum\_dictated
- Addendum\_Transcribing
- Addendum\_Transcribed
- Addendum\_Approved
- Read
- Config..

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